



SPA – Società Prodotti Antibiotici S.p.A.

Via Biella, 8 – 20143 Milano – Tel. +390289139.1 - Fax +39028132983 - email: biospa@spafarma.com

Product: **BIOTINYLATED GOAT ANTI-RABBIT IgG +
ACETYL-AVIDIN:BIOTINYLATED PEROXIDASE COMPLEX**
+ Normal Goat Serum

Code: **GR11-6**

Lot:

Description: Kit for immunohistochemistry complete of:
- normal goat serum
- biotinylated affinity-purified goat anti-rabbit IgG (H+L)
- pre-formed complex of modified avidin and biotinylated peroxidase

Appearance: Stabilized buffered stock solutions, to be diluted before use.

Manufacturing date:

PRODUCT PROFILE

... ml **Normal Goat Serum** Lot:
Suggested dilution 1: 30

... ml **Biotinylated affinity-purified goat anti-rabbit IgG** Lot:
Suggested dilution 1: 250

- a) Purity: this product was isolated from antisera by immunoaffinity chromatography using antigens coupled to agarose beads.
- b) Antibody specificity: the antibody used for biotinylation reacts with the heavy chains on rabbit IgG and with light chains common to most rabbit immunoglobulins. No antibody was detected against non-immunoglobulin serum proteins. This antibody has been tested to ensure minimal cross-reaction with human serum proteins, but the antibody may cross-react with immunoglobulins from other species.

..... ml **Modified avidin and biotinylated peroxidase complex** Lot:
Suggested dilution 1: 250

- a) Acetyl - avidin
Activity (biotin binding): 50 ± 15 units/ml solution.
Unit definition: one unit of avidin will bind 1 μ g of biotin.
- b) Peroxidase concentration: (spectrophot. absorbt. at 403 nm; an extinction coefficient of 18.6 has been used): 2500 μ g/ml solution.

Recommended storage conditions: 2-8 °C Stability: 12 months

FOR GENERAL LABORATORY USE